

PrP Antibody (R35863)

Catalog No.	Formulation	Size
R35863-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5621
Localization	Cytoplasmic
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PrP antibody is available for research use only.



Western blot testing of human brain lysate with PrP antibody at 0.3ug/ml. Expected molecular weight: 20~29kDa.

Description

Additional name(s) for this target protein: Prion protein, PRNP

Application Notes

Optimal dilution of the PrP antibody should be determined by the researcher.

Immunogen

Amino acids SDYEDRYREN were used as the immunogen for this PrP antibody.

Storage

Aliquot and store the PrP antibody at -20oC.