


PRRX1 Antibody (R35183)

Catalog No.	Formulation	Size
R35183-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5396
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PRRX1 antibody is available for research use only.



Western blot testing of human heart lysate with PRRX1 antibody at 1ug/ml. Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Paired related homeobox 1, PHOX1, PMX1

Application Notes

Optimal dilution of the PRRX1 antibody should be determined by the researcher.

Immunogen

Amino acids QADENVGEAGRS were used as the immunogen for this PRRX1 antibody.

Storage

Aliquot and store the PRRX1 antibody at -20oC.