


PSD95 Antibody / DLG4 (R33758)

Catalog No.	Formulation	Size
R33758-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1742
Localization	Cytoplasmic
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PSD95 antibody is available for research use only.



Western blot testing of rat brain lysate with PSD95 antibody at 0.3ug/ml. Molecular weight: routinely observed at 80~95 kDa.

Description

Additional name(s) for this target protein: Post-synaptic density protein 95; DLG4

Application Notes

Optimal dilution of the PSD95 antibody should be determined by the researcher.

Immunogen

Amino acids EEQARKAFDRATK were used as the immunogen for this PSD95 antibody.

Storage

Aliquot and store the PSD95 antibody at -20°C.