


PSMB4 Antibody (R34898)

Catalog No.	Formulation	Size
R34898-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Cow, Rabbit
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5692
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PSMB4 antibody is available for research use only.



Western blot testing of HEK293 lysate with PSMB4 antibody at 0.1ug/ml. Predicted molecular weight ~29 kDa, commonly observed at 26-29 kDa.



Western blot of mouse liver lysate with PSMB4 antibody at 0.3ug/ml. Predicted molecular weight ~29 kDa, commonly observed at 26-29 kDa.

Description

Additional name(s) for this target protein: Proteasome (prosome, macropain) subunit beta type 4, HN3, HsN3, PROS26, macropain beta chain

Application Notes

Optimal dilution of the PSMB4 antibody should be determined by the researcher.

Immunogen

Amino acids NRFQIATVTE were used as the immunogen for this PSMB4 antibody.

Storage

Aliquot and store the PSMB4 antibody at -20oC.