


## PSMB4 Antibody (R35302)

Catalog No.	Formulation	Size
R35302-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Cow, Rabbit
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5692
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PSMB4 antibody is available for research use only.



Western blot testing of HEK293 lysate with PSMB4 antibody at 0.3ug/ml. Predicted molecular weight ~29 kDa, commonly observed at 26-29 kDa.

### Description

Additional name(s) for this target protein: Proteasome beta 4 subunit

### Application Notes

Optimal dilution of the PSMB4 antibody should be determined by the researcher.

## Immunogen

Amino acids QIATVTEKGVE were used as the immunogen for this PSMB4 antibody.

## Storage

Aliquot and store the PSMB4 antibody at -20oC.