


PSMB7 Antibody (R35684)

Catalog No.	Formulation	Size
R35684-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5695
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PSMB7 antibody is available for research use only.



Western blot testing of HEK293 lysate with PSMB7 antibody at 0.3ug/ml. Predicted molecular weight: ~30kDa, observed here at ~40kDa.

Description

Additional name(s) for this target protein: Proteasome beta 7

Application Notes

Optimal dilution of the PSMB7 antibody should be determined by the researcher.

Immunogen

Amino acids RNAVLEADFAKR were used as the immunogen for this PSMB7 antibody.

Storage

Aliquot and store the PSMB7 antibody at -20oC.