


PSMB9 Antibody / Proteasome subunit beta type 9 (R35691)


Catalog No.	Formulation	Size
R35691-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Mouse, Dog, Pig, Cow, Rabbit
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5698
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PSMB9 antibody is available for research use only.



Western blot testing of human thymus lysate with PSMB9 antibody at 1ug/ml. Predicted molecular weight: ~23kDa.



Western blot of rat thymus lysate with PSMB9 antibody at 0.3ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: Proteasome beta 9

Application Notes

Optimal dilution of the PSMB9 antibody should be determined by the researcher.

Immunogen

Amino acids RNISYKYRED were used as the immunogen for this PSMB9 antibody.

Storage

Aliquot and store the PSMB9 antibody at -20oC.