


PTCD2 Antibody (R36207)

Catalog No.	Formulation	Size
R36207-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79810
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PTCD2 antibody is available for research use only.



Western blot testing of Daudi lysate with PTCD2 antibody at 1ug/ml. Predicted molecular weight: ~44kDa.

Description

Additional name(s) for this target protein: Pentatricopeptide repeat domain 2

Application Notes

Optimal dilution of the PTCD2 antibody should be determined by the researcher.

Immunogen

Amino acids TLHITGQVTDSL were used as the immunogen for this PTC2 antibody.

Storage

Aliquot and store the PTC2 antibody at -20°C.