


Ptgds Antibody (R36401)

Catalog No.	Formulation	Size
R36401-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	19215
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Ptgds antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with Ptgds antibody at 0.1 ug/ml.
Predicted molecular weight: ~21/27kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Prostaglandin D2 synthase

Application Notes

Optimal dilution of the Ptgds antibody should be determined by the researcher.

Immunogen

Amino acids YSRTQTLKDELKEK were used as the immunogen for this Ptgds antibody.

Storage

Aliquot and store the Ptgds antibody at -20oC.