


PTGIS Antibody / CYP8A1 (R36431)

Catalog No.	Formulation	Size
R36431-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5740
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PTGIS antibody is available for research use only.



Description

Additional name(s) for this target protein: Prostaglandin I2 synthase; Cytochrome P450, family 8, subfamily A, polypeptide 1

Application Notes

Optimal dilution of the PTGIS antibody should be determined by the researcher.

Immunogen

Amino acids RRAHRSKWLES were used as the immunogen for this PTGIS antibody.

Storage

Aliquot and store the PTGIS antibody at -20oC.