


## PTGIS Antibody / CYP8A1 (R36431)

Catalog No.	Formulation	Size
R36431-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Pig
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	5740
<b>Applications</b>	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This PTGIS antibody is available for research use only.



Western blot testing of human 1) lung and 2) ovary lysate with PTGIS antibody at 0.5ug/ml. Expected molecular weight: ~57 kDa.

### Description

Additional name(s) for this target protein: Prostaglandin I2 synthase; Cytochrome P450, family 8, subfamily A, polypeptide 1

### Application Notes

Optimal dilution of the PTGIS antibody should be determined by the researcher.

## **Immunogen**

Amino acids RRAHRSKWLES were used as the immunogen for this PTGIS antibody.

## **Storage**

Aliquot and store the PTGIS antibody at -20oC.