


PTGR1 Antibody / Prostaglandin reductase 1 (R36440)

Catalog No.	Formulation	Size
R36440-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22949
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PTGR1 antibody is available for research use only.



Western blot testing of human kidney lysate with PTGR1 antibody at 0.3ug/ml. Predicted molecular weight: ~36 kDa.

Description

Additional name(s) for this target protein: Prostaglandin reductase 1; PGR1; PRG-1; ZADH3

Application Notes

Optimal dilution of the PTGR1 antibody should be determined by the researcher.

Immunogen

Amino acids PTNSDFELKTAELP were used as the immunogen for this PTGR1 antibody.

Storage

Aliquot and store the PTGR1 antibody at -20oC.