


# PTGR2 Antibody / Prostaglandin reductase 2 (R36438)

Catalog No.	Formulation	Size
R36438-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	145482
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PTGR2 antibody is available for research use only.



#### **Description**

Additional name(s) for this target protein: Prostaglandin reductase 2; PRG-2; ZADH1; Zinc-binding alcohol dehydrogenase domain-containing protein 1

## **Application Notes**

Optimal dilution of the PTGR2 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids PDNINEGQVQVRT were used as the immunogen for this PTGR2 antibody. **Storage** Aliquot and store the PTGR2 antibody at -20oC.