


PTGR2 Antibody / Prostaglandin reductase 2 (R36438)

Catalog No.	Formulation	Size
R36438-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	145482
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PTGR2 antibody is available for research use only.



Western blot testing of pig testis lysate with PTGR2 antibody at 1ug/ml. Predicted molecular weight: ~38 kDa.

Description

Additional name(s) for this target protein: Prostaglandin reductase 2; PRG-2; ZADH1; Zinc-binding alcohol dehydrogenase domain-containing protein 1

Application Notes

Optimal dilution of the PTGR2 antibody should be determined by the researcher.

Immunogen

Amino acids PDNINEGQVQVRT were used as the immunogen for this PTGR2 antibody.

Storage

Aliquot and store the PTGR2 antibody at -20oC.