


PTP4A1 Antibody (R34384)

Catalog No.	Formulation	Size
R34384-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7803
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PTP4A1 antibody is available for research use only.



Western blot testing of A549 lysate with PTP4A1 antibody at 1ug/ml. Predicted/observed molecular weight: ~20kDa.

Description

Additional name(s) for this target protein: Protein tyrosine phosphatase type IVA

Application Notes

Optimal dilution of the PTP4A1 antibody should be determined by the researcher.

Immunogen

Amino acids KDSNGHRNNCCIQ were used as the immunogen for this PTP4A1 antibody.

Storage

Aliquot and store the PTP4A1 antibody at -20oC.