

PTPRC Antibody / CD45 [clone MSVA-045R] (V5980)

Catalog No.	Formulation	Size
V5980-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5980-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	MSVA-045R
UniProt	P08575
Localization	Cell membrane
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This PTPRC/CD45 antibody is available for research use only.

Application Notes

1. Optimal dilution of the PTPRC/CD45 antibody should be determined by the researcher.
2. This PTPRC/CD45 antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

A recombinant fragment (around amino acids 1-600) of human PTPRC protein (extracellular domain) was used as the immunogen for the PTPRC/CD45 antibody.

Storage

PTPRC/CD45 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

