


PUF60 Antibody (R36352)

Catalog No.	Formulation	Size
R36352-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22827
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PUF60 antibody is available for research use only.



Western blot testing of HEK293 lysate with PUF60 antibody at 0.1 ug/ml.
Predicted/observed molecular weight: ~60kDa.

Description

Additional name(s) for this target protein: Poly-U binding splicing factor 60KDa; SIAHBP1; FIR; fuse-binding protein-interacting repressor

Application Notes

Optimal dilution of the PUF60 antibody should be determined by the researcher.

Immunogen

Amino acids YDQERFDNSDLA were used as the immunogen for this PUF60 antibody.

Storage

Aliquot and store the PUF60 antibody at -20°C.