


PUM2 Antibody (R33599)

Catalog No.	Formulation	Size
R33599-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23369
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PUM2 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human PUM2-FLAG with PUM2 antibody (0.5ug/ml) in Lane A and probed with anti-FLAG (1/3000) in lane C. Mock-transfected HEK293 probed with PUM2 antibody (1ug/ml) in Lane B.

Description

Additional name(s) for this target protein: Pumilio homolog 2

Application Notes

Optimal dilution of the PUM2 antibody should be determined by the researcher.

Immunogen

Amino acids DQNRDLKQGDDDDSK were used as the immunogen for this PUM2 antibody.

Storage

Aliquot and store the PUM2 antibody at -20oC.