


Puma Antibody (F40032)


Catalog No.	Formulation	Size
F40032-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40032-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9BXH1
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This Puma antibody is available for research use only.



Puma antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human cervix tissue.



Puma antibody western blot analysis in A549 lysate. Predicted molecular weight: 21kDa.

Description

Puma is an essential mediator of p53-dependent and p53-independent apoptosis.

Application Notes

Titration of the Puma antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 109-137 from the alpha isoform was used as the immunogen for this Puma antibody.

Storage

Aliquot the Puma antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.