

PVALB Antibody / Calcium binding protein parvalbumin [clone PVALB/12794R] (V5983)

Catalog No.	Formulation	Size
V5983-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5983-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V5983SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	PVALB/12794R
UniProt	P20472
Localization	Cell junction, Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml Western Blot : 2-4ug/ml
Limitations	This PVALB/Calcium binding protein parvalbumin antibody is available for research use only.

Application Notes

- Optimal dilution of the PVALB/Calcium binding protein parvalbumin antibody should be determined by the researcher.
- This PVALB/Calcium binding protein parvalbumin antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

Recombinant full-length human PVALB protein was used as the immunogen for the PVALB/Calcium binding protein parvalbumin antibody.

Storage

PVALB/Calcium binding protein parvalbumin antibody with sodium azide - store at 2 to 8°C; antibody without sodium

azide - store at -20 to -80oC.