


## PYROXD1 Antibody (R36135)

Catalog No.	Formulation	Size
R36135-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79912
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PYROXD1 antibody is available for research use only.



Western blot testing of human colorectal cancer lysate with PYROXD1 antibody at 1ug/ml. Predicted/observed molecular weight ~56kDa. The ~48kDa band is likely a breakdown product.

## Description

Additional name(s) for this target protein: Pyridine nucleotide-disulphide oxidoreductase domain 1; Pyridine nucleotide-disulfide oxidoreductase domain-containing protein 1

## Application Notes

Optimal dilution of the PYROXD1 antibody should be determined by the researcher.

## Immunogen

Amino acids TCAEQLATHFPSED were used as the immunogen for this PYROXD1 antibody.

## Storage

Aliquot and store the PYROXD1 antibody at -20oC.