


QKI Antibody (R34017)

| Catalog No. | Formulation | Size |
|--------------|--|--------|
| R34017-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$ | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|----------------------|---|
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Cow, Dog, Pig, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 9444 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This QKI antibody is available for research use only. |



Description

Additional name(s) for this target protein: Quaking

Application Notes

Optimal dilution of the QKI antibody should be determined by the researcher.

Immunogen

Amino acids ETKEKPKPTPDY were used as the immunogen for this QKI antibody. **Storage** Aliquot and store the QKI antibody at -20oC.