


QKI Antibody (R34017)

Catalog No.	Formulation	Size
R34017-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9444
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This QKI antibody is available for research use only.



Western blot testing of mouse brain lysate with QKI antibody at 0.3ug/ml.
Predicted/observed molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Quaking

Application Notes

Optimal dilution of the QKI antibody should be determined by the researcher.

Immunogen

Amino acids ETKEKPKPTPDY were used as the immunogen for this QKI antibody.

Storage

Aliquot and store the QKI antibody at -20oC.