


## Rab17 Antibody (R35645)

Catalog No.	Formulation	Size
R35645-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	19329
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Rab17 antibody is available for research use only.



Western blot testing of mouse kidney lysate with Rab17 antibody at 0.2ug/ml.  
Predicted/observed molecular weight: ~23kDa.

## Description

Additional name(s) for this target protein: Ras related protein Rab17

## Application Notes

Optimal dilution of the Rab17 antibody should be determined by the researcher.

## **Immunogen**

Amino acids RKDSFHKAQQWLED were used as the immunogen for this Rab17 antibody.

## **Storage**

Aliquot and store the Rab17 antibody at -20oC.