


RAB2 Antibody (R34155)

Catalog No.	Formulation	Size
R34155-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5862
Applications	Western Blot : 1-3ug/ml
Limitations	This RAB2 antibody is available for research use only.



Western blot testing of HeLa lysate with RAB2 antibody at 1ug/ml. Isoforms of 24kDa and 21kDa can be detected.

Description

Additional name(s) for this target protein: GTP binding protein RAB2A; RAB2A

Application Notes

Optimal dilution of the RAB2 antibody should be determined by the researcher.

Immunogen

Amino acids HAGNQGGQQAGGGCC were used as the immunogen for this RAB2 antibody.

Storage

Aliquot and store the RAB2 antibody at -20oC.