


RAB23 Antibody (R35081)

Catalog No.	Formulation	Size
R35081-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51715
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RAB23 antibody is available for research use only.



Western blot testing of mouse eye lysate with RAB23 antibody at 0.1ug/ml. Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Ras-related protein Rab-23

Application Notes

Optimal dilution of the RAB23 antibody should be determined by the researcher.

Immunogen

Amino acids PNKQRTKKNRNP were used as the immunogen for this RAB23 antibody.

Storage

Aliquot and store the RAB23 antibody at -20oC.