

RAB27B Antibody (F42428)


Catalog No.	Formulation	Size
F42428-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42428-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O00194
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This RAB27B antibody is available for research use only.

95
72
55
36
28
17

RAB27B antibody western blot analysis in MDA-MB435 lysate.



RAB27B antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human stomach tissue.

Description

Members of the Rab protein family, including RAB27B, are prenylated, membrane-bound proteins involved in vesicular fusion and trafficking.[OMIM]

Application Notes

Titration of the RAB27B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 121-150 from the human protein was used as the immunogen for this RAB27B antibody.

Storage

Aliquot the RAB27B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.