


RAB5B Antibody (F46636)

Catalog No.	Formulation	Size
F46636-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F46636-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Chicken
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P61020
Applications	Western Blot : 1:1000
Limitations	This RAB5B antibody is available for research use only.



72 55 - 36 28 - ◀	RAB5B antibody western blot analysis in mouse liver tissue lysate
17	
10	

Description

Involved in protein transport and likely involved in vesicular traffic (By similarity). [UniProt]

Application Notes

Titration of the RAB5B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 85-112 from the human protein was used as the immunogen for this RAB5B antibody.

Storage

Aliquot the RAB5B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.