


RAB5C Antibody (F42976)


Catalog No.	Formulation	Size
F42976-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42976-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P51148
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This RAB5C antibody is available for research use only.



RAB5C antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human testis tissue.



RAB5C antibody western blot analysis in Jurkat lysate.

Description

Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment.

Application Notes

Titration of the RAB5C antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 178-207 from the human protein was used as the immunogen for this RAB5C antibody.

Storage

Aliquot the RAB5C antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.