

RACGAP1 Antibody (Biotin Conjugate) (R34126BTN)

Catalog No.	Formulation	Size
R34126BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29127
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This RACGAP1 antibody is available for research use only.



Western blot testing of human Jurkat lysate with biotinylated RACGAP1 antibody at 1ug/ml. Predicted molecular weight: ~71kDa.

Description

Additional name(s) for this target protein: MgcRacGAP, Rac GTPase activating protein 1

Application Notes

Optimal dilution of the RACGAP1 antibody should be determined by the researcher.

Immunogen

Amino acids GRQGNFFASPMLK were used as the immunogen for this RACGAP1 antibody.

Storage

Aliquot and store the RACGAP1 antibody at -20oC.