


RAD51AP1 Antibody (R33169)

Catalog No.	Formulation	Size
R33169-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10635
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RAD51AP1 antibody is available for research use only.



IHC testing of FFPE human testis with RAD51AP1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: RAD51-interacting protein

Application Notes

Optimal dilution of the RAD51AP1 antibody should be determined by the researcher.

Immunogen

Amino acids ARVKPLHPNATST were used as the immunogen for this RAD51AP1 antibody.

Storage

Aliquot and store the RAD51AP1 antibody at -20oC.