


RAD51C Antibody (R35618)

Catalog No.	Formulation	Size
R35618-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5889
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RAD51C antibody is available for research use only.



Western blot testing of HeLa lysate with RAD51C antibody at 2ug/ml. Predicted molecular weight: ~42kDa.

Description

Additional name(s) for this target protein: RAD51 homolog C, RAD51 homolog 3

Application Notes

Optimal dilution of the RAD51C antibody should be determined by the researcher.

Immunogen

Amino acids RGKTFRFEMQRDL were used as the immunogen for this RAD51C antibody.

Storage

Aliquot and store the RAD51C antibody at -20oC.