


Rad51D Antibody [clone 28R25] (RQ8288)

Catalog No.	Formulation	Size
RQ8288	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28R25
Purity	Affinity chromatography
UniProt	O75771
Applications	Western Blot : 1:500-1:2000
Limitations	This Rad51D antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with Rad51D antibody. Predicted molecular weight ~35 kDa.

Description

DNA repair protein RAD51 homolog 4 (Rad51D) is involved in the homologous recombination repair (HRR) pathway of double-stranded DNA breaks arising during DNA replication or induced by DNA-damaging agents. [UniProt]

Application Notes

Optimal dilution of the Rad51D antibody should be determined by the researcher.

Immunogen

A synthetic peptide derived from human Rad51D was used as the immunogen for the Rad51D antibody.

Storage

Store the Rad51D antibody at -20°C.