


## Ramp1 Antibody (R33685)

Catalog No.	Formulation	Size
R33685-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10267
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Ramp1 antibody is available for research use only.



Western blot testing of rat brain lysate with Ramp1 antibody at 0.3ug/ml. Routinely observed at ~18kDa and 40kDa.

## Description

Additional name(s) for this target protein: Receptor (calcitonin) activity modifying protein 1

## Application Notes

Optimal dilution of the Ramp1 antibody should be determined by the researcher.

1. The Ramp1 antibody is specific for isoform 1 of the protein.

## Immunogen

Amino acids DYQSKRTEGIV were used as the immunogen for this Ramp1 antibody.

## Storage

Aliquot and store the Ramp1 antibody at -20oC.