


RanBP1 Antibody (R34744)


Catalog No.	Formulation	Size
R34744-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5902
Localization	Cytoplasmic
Applications	Western Blot : 0.05-0.1ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RanBP1 antibody is available for research use only.



Western blot testing of HeLa lysate with RanBP1 antibody at 0.05ug/ml. Expected molecular weight: 23-26 kDa.



IHC staining of FFPE human kidney with RanBP1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: RAN binding protein 1

Application Notes

Optimal dilution of the RanBP1 antibody should be determined by the researcher.

Immunogen

Amino acids LSVKEETKEDAEKQ were used as the immunogen for this RanBP1 antibody.

Storage

Aliquot and store the RanBP1 antibody at -20oC.