


Rapsyn Antibody (R36393)

Catalog No.	Formulation	Size
R36393-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5913
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Rapsyn antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with Rapsyn antibody at 0.5ug/ml.
Predicted molecular weight: ~40kDa, routinely observed at 40-43kDa.

Description

Additional name(s) for this target protein: 43 kD Receptor-associated protein of the synapse; RAPSN

Application Notes

Optimal dilution of the Rapsyn antibody should be determined by the researcher.

Immunogen

Amino acids YQSNQTEKALQVWTK were used as the immunogen for this Rapsyn antibody.

Storage

Aliquot and store the Rapsyn antibody at -20oC.