


Rapsyn Antibody (R36393)

Catalog No.	Formulation	Size
R36393-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5913
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Rapsyn antibody is available for research use only.



Description

Additional name(s) for this target protein: 43 kD Receptor-associated protein of the synapse; RAPSN

Application Notes

Optimal dilution of the Rapsyn antibody should be determined by the researcher.

Immunogen

Amino acids YQSNQTEKALQVWTK were used as the immunogen for this Rapsyn antibody. **Storage** Aliquot and store the Rapsyn antibody at -20oC.