


RASSF2 Antibody / Ras association domain-containing protein 2 [clone 28R94] (RQ8517)

Catalog No.	Formulation	Size
RQ8517	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28R94
Purity	Affinity chromatography
UniProt	P50749
Applications	Western Blot : 1:500-1:2000
Limitations	This RASSF2 antibody is available for research use only.



Western blot testing of human Jurkat cell lysate with RASSF2 antibody. Predicted molecular weight ~38 kDa.

Description

Potential tumor suppressor. Acts as a KRAS-specific effector protein. May promote apoptosis and cell cycle arrest. Stabilizes STK3/MST2 by protecting it from proteasomal degradation. [UniProt]

Application Notes

Optimal dilution of the RASSF2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human RASSF2 protein was used as the immunogen for the RASSF2 antibody.

Storage

Store the RASSF2 antibody at -20oC.