


RASSF3 Antibody (R33265)

Catalog No.	Formulation	Size
R33265-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	283349
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This RASSF3 antibody is available for research use only.



Western blot testing of heart lysate with RASSF3 antibody at 0.5ug/ml with [B] and without [A] immunizing/blocking peptide. Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Ras association (RalGDS/AF-6) domain family 3

Application Notes

Optimal dilution of the RASSF3 antibody should be determined by the researcher.

Immunogen

Amino acids CKPPQTSPNSGK were used as the immunogen for this RASSF3 antibody.

Storage

Aliquot and store the RASSF3 antibody at -20oC.