


# RASSF7 Antibody (R33349)

Catalog No.	Formulation	Size
R33349-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8045
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This RASSF7 antibody is available for research use only.



# **Description**

Additional name(s) for this target protein: Ras association (RalGDS/AF-6) domain family 7

## **Application Notes**

Optimal dilution of the RASSF7 antibody should be determined by the researcher.

### **Immunogen**

Amino acids CTDLRGLELRVQRN were used as the immunogen for this RASSF7 antibody.

# **Storage** Aliquot and store the RASSF7 antibody at -20oC.