


RASSF7 Antibody (R33349)

Catalog No.	Formulation	Size
R33349-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8045
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This RASSF7 antibody is available for research use only.



Western blot testing of Jurkat cell lysate with RASSF7 antibody at 1 ug/ml. Predicted molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Ras association (RalGDS/AF-6) domain family 7

Application Notes

Optimal dilution of the RASSF7 antibody should be determined by the researcher.

Immunogen

Amino acids CTDLRGLELRVQRN were used as the immunogen for this RASSF7 antibody.

Storage

Aliquot and store the RASSF7 antibody at -20oC.