


RBM3 Antibody (R33619)

Catalog No.	Formulation	Size
R33619-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5935
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This RBM3 antibody is available for research use only.



Western blot testing of HepG2 lysate with RBM3 antibody at 0.01ug/ml. The expected ~17kDa band and the additional ~48kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: RNA binding motif (RNP1, RRM) protein 3

Application Notes

Optimal dilution of the RBM3 antibody should be determined by the researcher.

Immunogen

Amino acids DRYSGGNYRDNYDN were used as the immunogen for this RBM3 antibody.

Storage

Aliquot and store the RBM3 antibody at -20oC.