


RBMXL2 Antibody / HNRNPGT / RNA-binding motif protein X-linked-like-2 [clone 29R48] (RQ8748)

Catalog No.	Formulation	Size
RQ8748	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	29R48
Purity	Affinity chromatography
UniProt	O75526
Applications	Western Blot : 1:500-1:2000 Immunoprecipitation : 1:50 per 500ug of lysate
Limitations	This RBMXL2 antibody is available for research use only.



Immunoprecipitation of RBMXL2 protein from 500ug of human HeLa whole cell lysate using RBMXL2 antibody at 1:50 (Western blot done at 1:1000).



Western blot testing of human HeLa cell lysate with RBMXL2 antibody. Predicted molecular weight ~43 kDa.

Description

Additional names: RNA-binding motif protein X-linked-like-2; RBMXL2; Testis-specific heterogeneous nuclear ribonucleoprotein G-T; hnRNP G-T; HNRNPGT

Application Notes

Optimal dilution of the RBMXL2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human RBMXL2 protein was used as the immunogen for the RBMXL2 Antibody.

Storage

Store the RBMXL2 antibody at -20oC.