


## RBMXL2 Antibody / HNRNPGT / RNA-binding motif protein X-linked-like-2 [clone 29R48] (RQ8748)

Catalog No.	Formulation	Size
RQ8748	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

<b>Availability</b>	1-3 days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	29R48
<b>Purity</b>	Affinity chromatography
<b>UniProt</b>	O75526
<b>Applications</b>	Western Blot : 1:500-1:2000 Immunoprecipitation : 1:50 per 500ug of lysate
<b>Limitations</b>	This RBMXL2 antibody is available for research use only.



Immunoprecipitation of RBMXL2 protein from 500ug of human HeLa whole cell lysate using RBMXL2 antibody at 1:50 (Western blot done at 1:1000).



Western blot testing of human HeLa cell lysate with RBMXL2 antibody. Predicted molecular weight ~43 kDa.

## Description

Additional names: RNA-binding motif protein X-linked-like-2; RBMXL2; Testis-specific heterogeneous nuclear ribonucleoprotein G-T; hnRNP G-T; HNRNPGT

## Application Notes

Optimal dilution of the RBMXL2 antibody should be determined by the researcher.

## Immunogen

A synthetic peptide specific to human RBMXL2 protein was used as the immunogen for the RBMXL2 Antibody.

## Storage

Store the RBMXL2 antibody at -20°C.