


RBP4 Antibody / Retinol Binding Protein 4 [clone ACDC-18] (RQ5008)

Catalog No.	Formulation	Size
RQ5008	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	ACDC-18
Purity	Affinity purified
UniProt	P02753
Applications	Western Blot : 1:1000-1:5000
Limitations	This RBP4 antibody is available for research use only.



Western blot testing of human HepG2 cell lysate with RBP4 antibody. Predicted molecular weight ~23 kDa.

Description

Retinol-binding protein that mediates retinol transport in blood plasma. [UniProt]

Application Notes

Optimal dilution of the RBP4 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human RBP4 was used as the immunogen for the RBP4 antibody.

Storage

Store the RBP4 antibody at -20oC.