

RBP4 Antibody / Retinol Binding Protein 4 (R35327)

Catalog No.	Formulation	Size
R35327-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5950
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RBP4 antibody is available for research use only.



Western blot testing of human skin lysate with RBP4 antibody at 0.3ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: Retinol binding protein 4

Application Notes

Optimal dilution of the RBP4 antibody should be determined by the researcher.

Immunogen

Amino acids QKGNDDHWIVDTDYD were used as the immunogen for this RBP4 antibody.

Storage

Aliquot and store the RBP4 antibody at -20oC.