


Recombinant ACTN2 Antibody / Alpha Actinin 2 / Sarcomeric Alpha Actinin [clone rACTN2/7361] (V4782)

Catalog No.	Formulation	Size
V4782-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4782-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4782SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant MOUSE MONOCLONAL

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Mouse Monoclonal
Isotype	Mouse IgG1, kappa
Clone Name	rACTN2/7361
Purity	Protein A/G affinity
UniProt	P35609
Localization	Cytoplasm
Applications	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Limitations	This recombinant ACTN2 antibody is available for research use only.



IHC staining of FFPE human heart muscle tissue with recombinant ACTN2 antibody (clone rACTN2/7361). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.

Description

Alpha actinins belong to the spectrin gene superfamily which represents a diverse group of cytoskeletal proteins, including the alpha and beta spectrins and dystrophins. Alpha actinin is an actin-binding protein with multiple roles in different cell types. In non-muscle cells, the cytoskeletal isoform is found along microfilament bundles and adherens-type

junctions, where it is involved in binding actin to the membrane. In contrast, skeletal, cardiac, and smooth muscle isoforms are localized to the Z-disc and analogous dense bodies, where they help anchor the myofibrillar actin filaments. This gene encodes a muscle-specific, alpha actinin isoform that is expressed in both skeletal and cardiac muscles. Several transcript variants encoding different isoforms have been found for this gene.

Application Notes

Optimal dilution of the recombinant ACTN2 antibody should be determined by the researcher.

Immunogen

A portin of amino acids 557-692 from the human protein was used as the immunogen for the recombinant ACTN2 antibody.

Storage

Aliquot the recombinant ACTN2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.