


Recombinant Aurora B Antibody / AURKB [clone AURKB/3121R] (V7522)

Catalog No.	Formulation	Size
V7522-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V7522-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V7522SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug
V7522IHC-7ML	Prediluted in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide; *For IHC use only*	7 ml


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	AURKB/3121R
Purity	Protein A affinity
UniProt	Q96GD4
Localization	Nuclear, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This recombinant Aurora B antibody is available for research use only.



IHC testing of FFPE human colon tissue with recombinant Aurora B antibody (clone AURKB/3121R). HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min followed by cooling at RT for 20 min.



SDS-PAGE analysis of purified, BSA-free recombinant Aurora B antibody (clone AURKB/3121R) as confirmation of integrity and purity.

Description

Recognizes a protein of 39kDa, which is identified as Aurora B. The serine/threonine protein kinase aurora B (Aurora B) is a chromosomal passenger protein critical for accurate chromosome segregation, cytokinesis, protein localization to the centromere and kinetochore, correct microtubule-kinetochore attachment, and regulation of the mitotic checkpoint. Aurora B forms a tight complex with inner centrosome protein and survivin. Inactivation of any of these proteins causes similar defects in chromosome segregation. A significant overexpression of Aurora B has been found in a variety of human tumors including non-small cell lung carcinoma, astrocytoma, seminoma and carcinomas of the colon, prostate, endometrium and thyroid. The expression level of Aurora B is associated with cell proliferation and prognosis in these tumors.

Application Notes

Optimal dilution of the recombinant Aurora B antibody should be determined by the researcher.

Immunogen

A portion of amino acids 89-251 was used as the immunogen for this recombinant Aurora B antibody.

Storage

Store the recombinant Aurora B antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).