


Recombinant Bax Antibody [clone BAX/8960R] (V4648)

Catalog No.	Formulation	Size
V4648-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4648-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4648SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	BAX/8960R
Purity	Protein A/G affinity
UniProt	Q07812
Localization	Cytoplasm
Applications	Western Blot : 1-2ug/ml
Limitations	This recombinant Bax antibody is available for research use only.



Western blot testing of human A549 cell lysate with recombinant Bax antibody (clone BAX/8960R). Predicted molecular weight ~21 kDa.



SDS-PAGE analysis of purified, BSA-free Bax antibody (clone BAX/8960R) as confirmation of integrity and purity.

Description

Recognizes a protein of 21kDa, identified as the Bax protein. This mAb is highly specific to Bax and shows no cross-reaction with Bcl-2 or Bcl-X protein. Bcl-2 blocks cell death following a variety of stimuli. Bax has extensive amino acid homology with Bcl-2 and it homodimerizes and forms heterodimers with Bcl-2. Overexpression of Bax accelerates apoptotic death induced by cytokine deprivation in an IL-3 dependent cell line, and Bax also counters the death repressor activity of Bcl-2.

Application Notes

Optimal dilution of the recombinant Bax antibody should be determined by the researcher.

Immunogen

A portion of amino acids 1-100 from human BAX protein was used as the immunogen for the recombinant Bax antibody.

Storage

Aliquot the recombinant Bax antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.