


Recombinant Catenin gamma Antibody / Plakoglobin [clone CTNG/2155R] (V7899)

Catalog No.	Formulation	Size
V7899-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V7899-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V7899SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	CTNG/2155R
Purity	Protein A affinity chromatography
UniProt	P14923
Applications	ELISA (order BSA-free Format For Coating) : Flow Cytometry : 1-2ug/10 ⁶ cells in 0.1ml Immunofluorescence : 1-2ug/ml Western Blot : 1-2ug/ml
Limitations	This recombinant Catenin gamma antibody is available for research use only.



SDS-PAGE analysis of purified, BSA-free recombinant Catenin gamma antibody (clone CTNG/2155R) as confirmation of integrity and purity.

Description

It recognizes a protein of 80-87kDa, identified as gamma-catenin. The catenins (-T-catenin, is located on chromosome 10, and mutations in this gene show a strong correlation to late-onset Alzheimer's disease (LOAD) as well as to dilated cardiomyopathy.

Application Notes

Optimal dilution of the recombinant Catenin gamma antibody should be determined by the researcher.

Immunogen

Recombinant human full length protein was used as the immunogen for the recombinant Catenin gamma antibody.

Storage

Store the recombinant Catenin gamma antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).